

TRAVERSA**BRODA[®]**

Introducing the Traversa™ Transport Wheelchair – a cutting-edge innovation in non-emergency medical transportation (NEMT). Meticulously designed and crash-tested to meet the highest standards of comfort and safety, the Traversa is the first wheelchair solution made especially for NEMT. Backed by Broda's 40+ years of experience in wheelchair comfort, the Traversa is redefining the way passengers prefer to travel.



OPERATOR BENEFITS


- Replaces up to 90% of your stretcher trips
- Increased comfort during transit
- Dignified travel experience
- Lower risk of injury
- Simple adjustments for a customized wheelchair fit
- Easy transfers in and out of the wheelchair
- Faster service to and from the destination

**PASSENGER BENEFITS**

- Saves money on labour and fleet management
- Reduces injury risk, with no manual lifting required
- Quicker pickups for time-sensitive trips
- Maneuvers easily in tight spaces
- Crash-tested and fully WC19 compliant
- Integrates seamlessly with most vehicles
- Rapid return on investment



**SAVE TIME & TIME****ENHANCE PASSENGER
COMFORT****REDUCE RISK****CONTACT US**


WA  (08) 9456 1661
 saleswa@hospequip.com.au

VIC  (03) 9580 1055
 vicadmin@hospequip.com.au

 hospequip.com.au

SA  (08) 7225 8720
 sareception@hospequip.com.au

QLD  (07) 3448 9243
 qldadmin@hospequip.com.au

NSW  1300 467 737
 nswadmin@hospequip.com.au



CRASH-TESTED



REDUCES WORKLOAD



IMPROVES SAFETY



ADJUSTABLE PARTS





THE TRAVERSA TRANSPORT WHEELCHAIR CAN HANDLE UP TO 90% OF STRETCHER TRIPS.



	Traversa™	Stretcher
Made for NEMT	✓	✗
Comfortable for Passengers	✓	✗
Adjustable for Many Passenger Sizes	✓	✗
Can be used by a single operator	✓	✗
Adjustable Components	✓	✗
Improves Safety	✓	✗
Affordable	✓	✗
High Return on Investment	✓	✗





Overall Product Dimensions	Degrees of Tilt	Degrees of Recline	SWL	Adjustable Seat Widths	Seat Depth	Back Height	Food Pump Adjustable Seat Surface Height	Warranty
701mm W x 1320mm H	Up to 20°	Up to 90°	204Kg / 450lbs	51 - 61cm 20 - 24"	47cm 18.5"	82cm 32"	57cm - 81cm 22.5 - 32"	5 Years on Frame


CONTACT US



WA  (08) 9456 1661
 saleswa@hospequip.com.au

VIC  (03) 9580 1055
 vicadmin@hospequip.com.au

 hospequip.com.au

SA  (08) 7225 8720
 sareception@hospequip.com.au

QLD  (07) 3448 9243
 qldadmin@hospequip.com.au

NSW  1300 467 737
 nswadmin@hospequip.com.au